## Ballancer Aesthetic Line



Ballancer® systems have introduced a revolution in aesthetic compression treatment. Mego Afek has over 30 years of experience treating patients with its famous medical line of Lympha Press® products. All this experience is implemented in the clinically and scientifically proven **Ballancer**<sup>®</sup> systems.

Ballancer® treatments stimulate dermal and subdermal circulation. The increased blood flow and decongestive directional massage work together to tone and invigorate the skin, speed removal of fatiguing waste products, slim contours and create a feeling of relaxation and well-being.

Mego Afek AC Ltd. quality system complies with: Ballancer® 505 is FDA cleared and meets CE standards.





# Ballancer® 505V2

### Aesthetic Compression Treatment



Phiebo Press

## How the Ballancer<sup>®</sup> Works

Ballancer® uses a special inflating compression garment to apply a gentle or firm massage. Each Ballancer® garment contains up to 24 individual air chambers that overlap to apply fluent compression strokes that target each centimeter of the treated area. This massage motion, with its special direction of flow, increases circulation and gently but thoroughly decongests the problem areas of hips, thighs, buttocks and upper arms.

The Ballancer® action accelerates the removal of waste products from fat metabolism, reducing the appearance of cellulite ("orange peel skin"), improving skin tone, and reducing volume in areas affected by excess fluid retention. Its gentle, rhythmic action promotes relaxation.

### How does Ballancer<sup>®</sup> Affect Your Body?

#### **CIRCULATORY SYSTEM:**

Ballancer® peristaltic massage enhances blood circulation in the skin and the deep layers under it. This brings fresher oxygen to skin cells, and assists in speeding the elimination of toxins.

#### **DIGESTIVE SYSTEM:**

Ballancer<sup>®</sup> compresses and massages the abdomen, relaxing the abdominal muscles. The massage action is a natural bowel stimulant and can be especially helpful for clients who suffer from constipation.

#### LYMPHATIC SYSTEM:

Ballancer® increases the lymphatic flow. The lymphatic system is responsible for removal of waste products from the areas under the skin, and Ballancer® massage speeds their elimination.

#### **MUSCLES:**

Ballancer® peristaltic massage helps remove toxins from exercise and physical activity, as well as soothing, relaxing and decongesting the muscles, easing aches and pains.

#### SKIN:

**Ballancer**<sup>®</sup> massage increases the flow of freshly oxygenated blood to skin layers. Skin will look firm and toned, and will feel smoother and softer.

#### **Ballancer® Treatment Benefits**

- Treats and reduces the appearance of cellulite
- Reduces volume in areas of excess water retention, visibly slimming and toning contours
- Firms and tones the skin, for fresh, young look
- Clients complete their session feeling relaxed and rejuvenated

# Ballancer® 505V2


- Ballancer® 505V2, with its unique features and simple use, is ideal for beautysalons, studios and spas.
- Ballancer® 505V2 provides these treatment modes: Pre-treatment, Ballancer<sup>®</sup>, Wave
- The Pretreatment cycle is a primer for the main treatment. It decongests the trunk first, to maximize the results of the t session that follows. As the sequence progresses, it begins compression moving lower and lower down the limb.
- Ballancer® applies two treatment modes, according to sele
  - Wave cycle long, deep, and slow peristaltic massage
  - Ballancer® cycle rapid, superficial peristaltic caresses

### **Unique Garments**

Our innovative, unique and patent protected garments are world -famous for their high quality, optimal treatment and comfort.

- · Garments that treat entire upper or lower body
- Patent protected, fully overlapping chambers unique system for smooth and effective treatment, without gaps.
- More compression cells 24 overlapping compression chambers for efficient, fluent, comfortable massage
- Ballancer® Pants can accommodate most sizes due a special design !
- · Sturdy and durable fabrics that are easy to maintain and keep clean







	Ballancer <sup>®</sup> 505	
	Technical specifications:	
	Model	1201
treatment the	Cycle	Pretreatment, Ballancer® Wave
	Pressure Range	20 - 80 mmHg
ection:	Size	19x 31x 38 cm 7.5"x 12.2"X 14.9"
	Weight	6.7 Kg ( 14.7 lbs)
5	Electric Voltage	115V - 230V 50/60 Hz